

# App Guide Hip Extension

SageMotion  
Wearable Biofeedback System





## **Table of Contents**

Components	<i>(page 1)</i>
Wirelessly Connect to Computer or Cellphone	<i>(page 2)</i>
Hip Extension App	<i>(page 3)</i>

# Components



Hub



Nodes (8x)



Battery



Node Straps: *Medium (8x), Short (4x), Long (2x)*



Cable A (10x)

-Connect Hub to Battery  
-Charge Nodes & Battery



Cable B (*optional use*)

-Connect Hub to Computer




Node Charging Station

# Wirelessly Connect to Computer or Cellphone

## 1) Connect Cable A to Battery and to Hub




## 2) On Computer/Cellphone, Connect to Wi-Fi: "SageMotion"



*Note 1: Need to wait for up to 1 minute for "SageMotion" to appear in Wi-Fi list. If it doesn't appear, try turning the Wi-Fi off and then on again on the computer/cellphone.*

*Note 2: Hub is connected after clicking "Connect" even if in Windows it shows "Connecting" or "No internet, open".*

## 3) On Computer/Cellphone, in Chrome Address Bar, Go To <http://192.168.12.1>



**[Note] If Computer Doesn't Have Wi-Fi:** plug in Cable B to the Hub and to the ethernet port of your computer, then in chrome address bar, go to **<http://192.168.137.1>**

# Hip Extension App

*The purpose of the Hip extension App is to record, analyze, and provide feedback for Hip Extension angle while subjects perform walking or other daily activities.*

## 1) Turn on 4 Nodes




*Slide switch toward middle to turn node on*



*Green light will blink after the node is on and running*

## 2) Select “Hip extension” App



## 3) Click “Search”

**Node List**





# Hip Extension App (cont.)

4) Configure 2 Sensor Nodes and 2 Feedback Nodes as Shown Below:

Node List

Search

Connect

Type	Position	MAC	
feedback ▾	feedback_max ▾	88:6B:0F:E1:D8:96	
sensor ▾	pelvis ▾	88:6B:0F:E1:D8:9E	
sensor ▾	thigh ▾	88:6B:0F:E1:D8:A2	
feedback ▾	feedback_min ▾	88:6B:0F:E1:D8:9F	

5) Click “Connect”

Node List

Search

Connect



6) “Ready to collect data” Will Appear after Node Connection is Complete

Hip Extension ▾


Start

Stop

✓ Ready to collect data


# Hip Extension App (cont.)

## 7) Thread Straps through Nodes and Attach at Locations Shown Below:



**[Note]** thigh (either left or right lateral, switch pointing up), pelvis (back center, switch pointing up). feedback\_min and feedback\_max nodes can be placed at any location

### How to Thread Straps



## 8) Click “Blink” for each Node to Confirm Correct Locations (red LED for given node blinks 3 times on click)

Type	Position	MAC		
sensor	thigh	88:6B:0F:E1:D8:A2		<input type="button" value="Blink"/>
sensor	pelvis	88:6B:0F:E1:D8:9E		<input type="button" value="Blink"/>
feedback	feedback_min	88:6B:0F:E1:D8:9F		<input type="button" value="Blink"/>
feedback	feedback_max	88:6B:0F:E1:D8:96		<input type="button" value="Blink"/>



# Hip Extension App (cont.)

9) In App Configuration, Enter Settings (Example Below)

## App Configuration

Trial Name	<input type="text" value="Hip_angle_1"/>
<b>Sensor Placement</b>	
Which Leg?	<input type="text" value="Left Leg"/> ▼
<b>Feedback Setting</b>	
Feedback On?	<input type="text" value="Yes"/> ▼
Min Threshold Angle	<input type="text" value="-20"/>
Max Threshold Angle	<input type="text" value="20"/>
<b>Save Options</b>	
Save Mode	<input type="text" value="xlsx"/> ▼

# Hip Extension App (cont.)

10) Click “Start” to Start Running the App





11) After the Trial is Finished, Click “Stop”





12) After Clicking “Stop”, a File from that Trial will Appear under Download Data. Click the File (e.g. Hip\_angle\_2) to Download it to the Computer or Phone.

## Data Management

 Download Selected

 Delete Selected

<input type="checkbox"/>	Name	Date▲	Duration	App	Type	Size	Rename	Delete
<input checked="" type="checkbox"/>	<u>Hip_angle_2</u>	2021-10-14-19-43-18	0:00:16	Hip Extension	.xlsx	1.8 MB		

# Hip Extension App (*cont.*)

## Description of Data in Downloaded File

**time** (sec): time since trial start

**Hip\_ext** (deg): Hip Flexion/Extension angle for the leg

**Feedback\_min**: Minimum threshold angle without feedback

**Feedback\_max**: Maximum threshold angle without feedback

**SensorIndex\_1/2**: index of raw sensor data

**AccelX/Y/Z\_1/2** (m/s<sup>2</sup>): raw acceleration data

**GyroX/Y/Z\_1/2** (deg/s): raw gyroscope data

**MagX/Y/Z\_1/2** (μT): raw magnetometer data

**Quat1/2/3/4\_1/2**: quaternion data

**Sampletime\_1/2**: timestamp of each sensor

**Package\_1/2**: package number of each sensor